

SUPER GRIP DISPOSABLE GLOVES FIT, DEXTERITY, TRACTION GRIP, AMBIDEXTROUS


- Embossed diamond pattern on the outside of the glove for ergonomic traction and exceptional wet grip.
- Resistant to Fentanyl
- Tested for use with Chemotherapy drugs
- Medical Grade
- Chemical resistant 100% 8 mil nitrile
- Beaded cuff for secure fit around the arm
- ◆ CE level 5 dexterity for superior fit, feel, and comfort
- Powder Free
- 9.5" Length
- Ambidextrous fit for less waste
- Industries: Industrial, Construction, Automotive, Agriculture



Tec	hnic	al D	ata

Material	100 % Nitrile	
Standards	ASTM D6978-05 - Resistance to Permeation of Permeation by Chemotherapy Drugs and Fentanyl ISO 374-1:2016 Type B - Chemical Resistance EN 374-2:2014 - Penetration Resistance EN 16523-1:2015 - Permeation Resistance EN 374-4:2013 - Resistance to Degradation by Chemicals EN ISO 374-5:2016 - Micro-organism Protection EN 455-1:2000 - Medical Glove for Single Use FDA CFR 21 177.2600 - Suitable for Food Contact	
Cuff Closure	Beaded Cuff	
Mil Thickness	8 Mil	
Colors	Black	
Sizes	S-XXL	
Packaging	100 gloves/box 10 boxes/carton	

Embossed diamond pattern provides superior grip





Ordering Information

Item	Size	UPC
3277DK/8	S	731791173900
3277DK/9	М	731791173917
3277DK/10	L	731791173924
3277DK/11	X1	731791173931
3277DK/12	X2	731791173948



















All safety equipment should be selected in conjunction with an approved safety program that assesses your specific hazards. The selection and use of specific products sold by Majestic is exclusively within the control of, and remains the sole responsibility of the customer. Many variables beyond the control of Majestic affect the results obtained in applying ese products. Inappropriate selection and/or use of safety products may cause serious injury.

SPECIFICATION SHEET 3277DK Reliability in Safety.

Majestic Glove 2510 W Casino Rd Everett, WA 98204

P: 800-367-4568 F: 425-407-1030 orders@majesticglove.com

majesticglove.com